

Product Datasheet

Guinea INS ELISA Kit (orb1197768)

Description	Guinea pig Insulin (INS) ELISA Kit
Reactivity	Guinea pig
Range	140-9000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01329
Sensitivity	64.3pg/mL
Entrez	100379579
Expiration Date	Please enquire.