

Product Datasheet

Mouse HSPA5 ELISA Kit (orb1197781)

Description	Mouse Heat Shock 70kDa Protein 5 (HSPA5) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P20029
Sensitivity	0.054ng/mL
Entrez	14828
Expiration Date	Please enquire.