

Product Datasheet

Mouse HSP90aA1 ELISA Kit (orb1197782)

Description	Mouse Heat Shock Protein 90kDa Alpha A1 (HSP90aA1) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q80Y52
Sensitivity	31pg/mL
Entrez	15519
Expiration Date	Please enquire.