

## Product Datasheet

### Mouse FOLR1 ELISA Kit (orb1197844)

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | Mouse Folate Receptor 1, Adult (FOLR1)<br>ELISA Kit           |
| <b>Reactivity</b>          | Mouse                                                         |
| <b>Range</b>               | 78.1-5000pg/mL                                                |
| <b>Tested Applications</b> | ELISA                                                         |
| <b>Note</b>                | For research use only                                         |
| <b>Application notes</b>   | sandwich                                                      |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>P35846</b>                                                 |
| <b>Sensitivity</b>         | 31pg/mL                                                       |
| <b>Entrez</b>              | <b>14275</b>                                                  |
| <b>Expiration Date</b>     | Please enquire.                                               |