

Product Datasheet

Mouse FGL1 ELISA Kit (orb1197846)

Description	Mouse Fibrinogen Like Protein 1 (FGL1) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q71KU9
Sensitivity	11.8pg/mL
Entrez	234199
Expiration Date	Please enquire.