

Product Datasheet

Mouse CRTAP ELISA Kit (orb1197918)

Description	Mouse Cartilage Associated Protein (CRTAP) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9CYD3
Sensitivity	0.056ng/mL
Entrez	56693
Expiration Date	Please enquire.