

Product Datasheet

Mouse COX3 ELISA Kit (orb1197931)

Description	Mouse Cytochrome C Oxidase Subunit III (COX3) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9MQW4
Sensitivity	0.119ng/mL
Entrez	3338893
Expiration Date	Please enquire.