

## Product Datasheet

### Mouse CHAF1B ELISA Kit (orb1197959)

|                            |                                                                 |
|----------------------------|-----------------------------------------------------------------|
| <b>Description</b>         | Mouse Chromatin Assembly Factor 1, Subunit B (CHAF1B) ELISA Kit |
| <b>Reactivity</b>          | Mouse                                                           |
| <b>Range</b>               | 0.156-10ng/mL                                                   |
| <b>Tested Applications</b> | ELISA                                                           |
| <b>Note</b>                | For research use only                                           |
| <b>Application notes</b>   | sandwich                                                        |
| <b>Sample Types</b>        | tissue homogenates, cell lysates or other biological fluids.    |
| <b>Uniprot ID</b>          | <b>Q9D0N7</b>                                                   |
| <b>Sensitivity</b>         | 0.052ng/mL                                                      |
| <b>Entrez</b>              | <b>110749</b>                                                   |
| <b>Expiration Date</b>     | Please enquire.                                                 |