

Product Datasheet

Mouse CDK8 ELISA Kit (orb1197964)

Description	Mouse Cyclin Dependent Kinase 8 (CDK8) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q8R3L8
Sensitivity	0.053ng/mL
Entrez	264064
Expiration Date	Please enquire.