

Product Datasheet

Canine CCK ELISA Kit (orb1197967)

Description	Canine Cholecystokinin (CCK) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9TS44
Sensitivity	4.64pg/mL
Entrez	609547
Expiration Date	Please enquire.