

Product Datasheet

Mouse CBX3 ELISA Kit (orb1197970)

Description	Mouse Chromobox Homolog 3 (CBX3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P23198
Sensitivity	0.046ng/mL
Entrez	12417
Expiration Date	Please enquire.