

Product Datasheet

Mouse ATP1b1 ELISA Kit (orb1198003)

Description	Mouse ATPase, Na ⁺ /K ⁺ Transporting Beta 1 Polypeptide (ATP1b1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P14094
Sensitivity	0.068ng/mL
Entrez	11931
Expiration Date	Please enquire.