

Product Datasheet

Mouse ATP1a2 ELISA Kit (orb1198004)

Description	Mouse ATPase, Na ⁺ /K ⁺ Transporting Alpha 2 Polypeptide (ATP1a2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q6PIE5
Sensitivity	0.042ng/mL
Entrez	98660
Expiration Date	Please enquire.