

Product Datasheet

Mouse AST2 ELISA Kit (orb1198008)

Description	Mouse Aspartate Aminotransferase 2 (AST2) ELISA Kit
Reactivity	Mouse
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P05202
Sensitivity	0.29ng/mL
Entrez	14719
Expiration Date	Please enquire.