

Product Datasheet

Mouse ARHGEF7 ELISA Kit (orb1198016)

Description	Mouse Rho Guanine Nucleotide Exchange Factor 7 (ARHGEF7) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9ES28
Sensitivity	34pg/mL
Entrez	54126
Expiration Date	Please enquire.