

Product Datasheet

Mouse ALDM ELISA Kit (orb1198032)

Description	Mouse Aldehyde Dehydrogenase, Mitochondrial (ALDM) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P47738
Sensitivity	0.12ng/mL
Entrez	11669
Expiration Date	Please enquire.