

Product Datasheet

Mouse ACP2 ELISA Kit (orb1198048)

Description	Mouse Acid Phosphatase 2, Lysosomal (ACP2) ELISA Kit
Reactivity	Mouse
Range	0.781-50U/L
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P24638
Sensitivity	0.37U/L
Entrez	11432
Expiration Date	Please enquire.