

Product Datasheet

Equine TGFb3 ELISA Kit (orb1198052)

Description	Equine Transforming Growth Factor Beta 3 (TGFb3) ELISA Kit
Reactivity	Equine
Range	125-8000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	F7DR45
Sensitivity	31pg/mL
Entrez	100051765
Expiration Date	Please enquire.