

Product Datasheet

Sheep TNFa ELISA Kit (orb1198054)

Description	Sheep Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Reactivity	Sheep
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P23383
Sensitivity	7.1pg/mL
Entrez	443540
Expiration Date	Please enquire.