

Product Datasheet

Monkey TGFb1 ELISA Kit (orb1198056)

Description	Monkey Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit
Reactivity	Monkey
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernates or other biological fluids
Sensitivity	13.4pg/mL
Entrez	574135
Expiration Date	Please enquire.