

Product Datasheet

Sheep IL6 ELISA Kit (orb1198060)

Description	Sheep Interleukin 6 (IL6) ELISA Kit
Reactivity	Sheep
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P29455
Sensitivity	6.8pg/mL
Entrez	443406
Expiration Date	Please enquire.