

Product Datasheet

Monkey IFN γ ELISA Kit (orb1198065)

Description	Monkey Interferon Gamma (IFN γ) ELISA Kit
Reactivity	Monkey
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P63309
Sensitivity	6.1pg/mL
Entrez	102128291
Expiration Date	Please enquire.