

Product Datasheet

Human TPP2 ELISA Kit (orb1198090)

Description	Human Tripeptidyl-Peptidase 2 (TPP2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P29144
Sensitivity	0.130ng/mL
Entrez	7174
Expiration Date	Please enquire.