

Product Datasheet

Rat Rat Vasopressin V1a Receptor (AVPR1A) ELISA Kit (orb1198098)

Description	Rat Vasopressin V1a Receptor (AVPR1A) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P30560
Sensitivity	0.130ng/mL
Entrez	25107
Expiration Date	Please enquire.