

Product Datasheet

Human TRPM1 ELISA Kit (orb1198188)

Description	Human Transient Receptor Potential Cation Channel Subfamily M, Member 1 (TRPM1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q7Z4N2
Sensitivity	0.051ng/mL
Entrez	4308
Expiration Date	Please enquire.