

Product Datasheet

Human TRAF6 ELISA Kit (orb1198195)

Description	Human TNF Receptor Associated Factor 6 (TRAF6) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9Y4K3
Sensitivity	0.053ng/mL
Entrez	7189
Expiration Date	Please enquire.