

Product Datasheet

Human TNNT1 ELISA Kit (orb1198215)

Description	Human Troponin T Type 1, Slow Skeletal (TNNT1) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P13805
Sensitivity	6.1pg/mL
Entrez	7138
Expiration Date	Please enquire.