

## Product Datasheet

### Human TNNC2 ELISA Kit (orb1198218)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | Human Troponin C Type 2, Fast (TNNC2)<br>ELISA Kit               |
| <b>Reactivity</b>          | Human                                                            |
| <b>Range</b>               | 31.2-2000pg/mL                                                   |
| <b>Tested Applications</b> | ELISA                                                            |
| <b>Note</b>                | For research use only                                            |
| <b>Application notes</b>   | sandwich                                                         |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates or<br>other biological fluids. |
| <b>Uniprot ID</b>          | <b>P02585</b>                                                    |
| <b>Sensitivity</b>         | 11.8pg/mL                                                        |
| <b>Entrez</b>              | <b>7125</b>                                                      |
| <b>Expiration Date</b>     | Please enquire.                                                  |