

## Product Datasheet

### Human TK1 ELISA Kit (orb1198231)

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | Human Thymidine Kinase 1, Soluble (TK1)<br>ELISA Kit          |
| <b>Reactivity</b>          | Human                                                         |
| <b>Range</b>               | 15.6-1000pg/mL                                                |
| <b>Tested Applications</b> | ELISA                                                         |
| <b>Note</b>                | For research use only                                         |
| <b>Application notes</b>   | sandwich                                                      |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>P04183</b>                                                 |
| <b>Sensitivity</b>         | 6.0pg/mL                                                      |
| <b>Entrez</b>              | <b>7083</b>                                                   |
| <b>Expiration Date</b>     | Please enquire.                                               |