

Product Datasheet

Human SUMO2 ELISA Kit (orb1198278)

Description	Human Small Ubiquitin Related Modifier Protein 2 (SUMO2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P61956
Sensitivity	0.101ng/mL
Entrez	6613
Expiration Date	Please enquire.