

Product Datasheet

Human STC1 ELISA Kit (orb1198291)

Description	Human Stanniocalcin 1 (STC1) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P52823
Sensitivity	27pg/mL
Entrez	6781
Expiration Date	Please enquire.