

Product Datasheet

Human SPTLC3 ELISA Kit (orb1198304)

Description	Human Serine Palmitoyltransferase, Long Chain Base Subunit 3 (SPTLC3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9NUV7
Sensitivity	0.058ng/mL
Entrez	55304
Expiration Date	Please enquire.