

Product Datasheet

Human SNAT1 ELISA Kit (orb1198323)

Description	Human Sodium Coupled Neutral Amino Acid Transporter 1 (SNAT1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9H2H9
Sensitivity	0.055ng/mL
Entrez	81539
Expiration Date	Please enquire.