

Product Datasheet

Human SIRT4 ELISA Kit (orb1198354)

Description	Human Sirtuin 4 (SIRT4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9Y6E7
Sensitivity	0.051ng/mL
Entrez	23409
Expiration Date	Please enquire.