

Product Datasheet

Human SIRT1 ELISA Kit (orb1198356)

Description	Human Sirtuin 1 (SIRT1) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q96EB6
Sensitivity	0.28ng/mL
Entrez	23411
Expiration Date	Please enquire.