

Product Datasheet

Human SAT1 ELISA Kit (orb1198403)

Description	Human Spermidine/Spermine N1-Acetyltransferase 1 (SAT1) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P21673
Sensitivity	0.29ng/mL
Entrez	6303
Expiration Date	Please enquire.