

Product Datasheet

Human S100A4 ELISA Kit (orb1198408)

Description	Human S100 Calcium Binding Protein A4 (S100A4) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P26447
Sensitivity	29pg/mL
Entrez	6275
Expiration Date	Please enquire.