

Product Datasheet

Human REG3g ELISA Kit (orb1198443)

Description	Human Regenerating Islet Derived Protein 3 Gamma (REG3g) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q6UW15
Sensitivity	14.5pg/mL
Entrez	130120
Expiration Date	Please enquire.