

Product Datasheet

Human QSOX1 ELISA Kit (orb1198463)

Description	Human Quiescin Q6 Sulfhydryl Oxidase 1 (QSOX1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	000391
Sensitivity	0.116ng/mL
Entrez	5768
Expiration Date	Please enquire.