

Product Datasheet

Human PLA2G2A ELISA Kit (orb1198557)

Description	Human Phospholipase A2, Group IIA (PLA2G2A) ELISA Kit
Reactivity	Human
Range	6.25-400pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P14555
Sensitivity	2.55pg/mL
Entrez	5320
Expiration Date	Please enquire.