

Product Datasheet

Human PAM ELISA Kit (orb1198604)

Description	Human Peptidylglycine Alpha Amidating Monooxygenase (PAM) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P19021
Sensitivity	13.5pg/mL
Entrez	5066
Expiration Date	Please enquire.