

Product Datasheet

Human PADI4 ELISA Kit (orb1198613)

Description	Human Peptidyl Arginine Deiminase Type IV (PADI4) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9UM07
Sensitivity	0.111ng/mL
Entrez	23569
Expiration Date	Please enquire.