

Product Datasheet

Human PADI2 ELISA Kit (orb1198614)

Description	Human Peptidyl Arginine Deiminase Type II (PADI2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9Y2J8
Sensitivity	0.112ng/mL
Entrez	11240
Expiration Date	Please enquire.