

Product Datasheet

Human OSGIN1 ELISA Kit (orb1198626)

Description	Human Oxidative Stress Induced Growth Inhibitor 1 (OSGIN1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9UJX0
Sensitivity	0.064ng/mL
Entrez	29948
Expiration Date	Please enquire.