

Product Datasheet

Human NEU ELISA Kit (orb1198697)

Description	Human Neuraminidase (NEU) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q99519
Sensitivity	0.126ng/mL
Entrez	4758
Expiration Date	Please enquire.