

Product Datasheet

Human MMRN1 ELISA Kit (orb1198762)

Description	Human Multimerin 1 (MMRN1) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q13201
Sensitivity	9.6pg/mL
Entrez	22915
Expiration Date	Please enquire.