

Product Datasheet

Human MHCB ELISA Kit (orb1198771)

Description	Human Major Histocompatibility Complex Class I B (MHCB) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P03989
Sensitivity	0.63ng/mL
Entrez	3106
Expiration Date	Please enquire.