

Product Datasheet

Human MDM2 ELISA Kit (orb1198783)

Description	Human Mdm2 p53 Binding Protein Homolog (MDM2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q00987
Sensitivity	0.064ng/mL
Entrez	4193
Expiration Date	Please enquire.