

Product Datasheet

Human MDH1 ELISA Kit (orb1198785)

Description	Human Malate Dehydrogenase 1 (MDH1) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P40925
Sensitivity	27pg/mL
Entrez	4190
Expiration Date	Please enquire.