

Product Datasheet

Human MCL1 ELISA Kit (orb1198787)

Description	Human Myeloid Cell Leukemia Sequence 1, Bcl2 Related (MCL1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q07820
Sensitivity	0.118ng/mL
Entrez	4170
Expiration Date	Please enquire.